

## Standardization of Potassium Permanganate



### General Aim

To determine the strength of potassium permanganate with a standard oxalic acid solution.

### Method

Determination of potassium permanganate by redox titration.

### Learning Objectives (ILOs)

- Define direct titration.
- Define the meaning of a standard solution.
- Differentiate between primary and secondary standards.
- Predict how a certain sample could be analyzed.
- Determine the strength of potassium permanganate with a standard oxalic acid solution.
- Perform titration with the control of the temperature.
- Calculate the molarity and the strength of a given standard based on a similar procedure.

### Theoretical Background/Context


- Quantitative analysis deals with the determination of the quantity of the substance to be analyzed. Methods of quantitative analysis may be classified into:
  - 1) **Gravimetric analysis:** It depends on isolating and weighing of the final product with known pure, stable and definite form.
  - 2) **Instrumental analysis:** It depends on measuring some physical properties which change quantitatively with changing concentration of sample.
  - 3) **Volumetric analysis (Titration):** It depends on measuring volume of standard solution (titrant) used for complete reaction with the sample.
- Titration is the capacity of the sample to combine with the suitable standard quantitatively through quantitative reaction.
- A quantitative reaction is the reaction that proceeds forward to produce stable product(s) such as weakly ionizable compounds, e.g. H<sub>2</sub>O, weak acids & base, sparingly soluble salts (precipitate), complex ion, etc. The types of quantitative Reactions can be:
  - 1) Neutralization Reactions
    - a. H<sub>2</sub>O formation
    - b. Displacement: Formation of weak acid or weak base
  - 2) Complexometric reactions
  - 3) Redox reactions (Electron transfer)
  - 4) Precipitometric reactions

## Theoretical Background/Context (Cont')

- Any sample is a solution of unknown concentration and a Standard is a solution of exactly known concentration. The requirements of titrimetric reactions are:

- The reaction must be simple and expressed by a chemical equation.
- A single reaction must occur between the sample and titrant.
- The reaction must be instantaneous (rapid).
- Suitable standard solutions must be available.
- The end point should be easily detected.

**Oxidation** is loss of electrons and increase in valency number, gain of oxygen or loss of hydrogen




**Reduction** is gain of electrons and reduction in valency number, loss of oxygen or gain of hydrogen



**Oxidizing agent or oxidant** is the substance that gain electrons



**Reducing agent or reductant** is the substance that donate (lose) electrons



## Principle of Work

- This experiment is a redox titration where potassium permanganate (KMnO<sub>4</sub>) is the titrant and oxalic acid is the analyte. Also, KMnO<sub>4</sub> is the oxidizing agent and oxalic acid is the reducing agent. The reaction between KMnO<sub>4</sub> and oxalic acid is carried out in an acidic medium because permanganate ion in the acidic medium is a very strong oxidizing agent. Acidity is introduced by adding dil. H<sub>2</sub>SO<sub>4</sub>.

- No indicators are used to determine the endpoint, because KMnO<sub>4</sub> is a self-indicator. Permanganate (MnO<sub>4</sub><sup>-</sup>) ion has a dark purple colour. In an acidic medium, MnO<sub>4</sub><sup>-</sup> is reduced to colorless manganous (Mn<sup>2+</sup>) ions. On reaching the end point, the addition of the last single drop of permanganate imparts a light purple colour to the solution.

